



Kit for semi-quantitative detection of Brettanomyces spp. in alcoholic beverages

Developed and patented by Grape srl

Description:

Self-brett® is an innovative device for semi-quantitative detection of Brettanomyces in alcoholic beverages. It is composed of a union between a species-specific culturedependent detection method and a simple ready-to-use technology (dip-slide) that permits the producer to perform a rapid self-controlled analysis.

Advantages:

- Autonomy: Self-brett® test can be carried out directly on the production site without involving laboratories or specialized people, while obtaining a clear result supported by multiple evidences.
- Clear result: the user can evaluate the presence of the spoilage microorganism in wine in three different ways:
 - Change of color (from blue to yellowish)
 - Develop of the typical 'Brett odor'
 - Observation of colonies on the agar medium







SMELL IF BRETT-ODORS DEVELOPED



OBSERVE IF YEAST COLONIES APPEARED IN THE AGAR

Cost: low compared to laboratory analysis.

Applications:

Evaluation for *Brettanomyces* presence in alcoholic beverages (wine, beer, cider)

Test Procedure:

- Unscrew and extract the slide from its cylindrical container, avoiding any contact with the agar surface
- Fill the container with the beverage to be analyzed until the designated spot on the label
- Shake gently and wait 60 sec
- Empty and drain
- Close and leave at room temperature for 7-10 days



Results interpretation:

If the tested sample is negative to all three conditions, then the result is negative (not contaminated).


+39 3935465396 E-mail self-brett@grapesrl.it

Grape Srl Gruppo Ricerche Avanzate Per l'Enologia Ampelion, Corso Enotria 2C ~ Alba (CN) P.iva: 03388640041





- If the medium color changes to yellowish, the typical 'Brett odor' develops and colonies are observed, the sample is contaminated.
- Based on the amount of colonies contamination entity can be estimated: **low, medium, high.**



Storage:

10-25°C away from light, until the expiration date on the label. Eliminate if signs of deterioration or contamination are evident.

Shelf-life:

6 months.

Quality control:

Each Self-brett® batch is tested with the following strains:

Strain Code	Species	Growth
DSM 20176	Lactobacillus hilgardii	NO
DSM 70001	Dekkera bruxellensis	YES
DSM 70451	Saccharomyces cerevisiae	NO
DSM 17891	Acetobacter oeni	NO

Tel +39 3935465396 E-mail self-brett@grapesrl.it

Grape Srl Gruppo Ricerche Avanzate Per l'Enologia Ampelion, Corso Enotria 2C ~ Alba (CN) P.iva: 03388640041

Frequently asked questions (FAQ):

- The color has changed but after 10 days there isn't any odor nor have colonies not developed
 in case of particularly acidic wine the agar medium can undergo a slight color variation at the moment of analysis.
- 2) I can't see the colonies, but the odor developed and the color has changed -> different Brettanomyces strains can show different morphology and colouring more or less evident, therefore a careful backlight observations advised.

Warnings and precautions:

The product is intended for professional use and must be used by properly trained operators only, in accordance with manufacturer's instruction. Self-brett® cannot be classified as hazardous product. For further information read the papers inside the package.

Grape srl will not be responsible anyhow for use non-compliancewith such indications. It is the responsibility of customers to act in compliance with the applicable health, environment and <u>safety</u> regulations and to take the necessary actions in relation to the storage, handling, sale and use of the products.

Disposal of waste:

Disposal of waste must be carried out according to national and local regulations in force. The product cannot be classified as a special hazardous waste according to the rules described in the decree (EC) 1357/2008.

The product is designed for professional use only

Attention: consult the document inside

The product has to be stored at indicated temperature

The product is disposable

Keep away from light